



Introduction

Name: ZL-1707 A series LCD display 0.01g high precision analytical balance



Tel: [+86-769-88055358](tel:+86-769-88055358)

Fax: [+86-769-88055357](tel:+86-769-88055357)

Website: www.zltester.com

E-mail: Zhongli@zltester.com

Address: [Room 722, Dongyuan Mansion, No.159 QiFeng Road, Dongguang city,P.R.C](#)



Specifications

- high precision analytical balance
- 1 Weighing Response time Adjustable
- 2 Counting And Weighing
- 3 Time & Date Setting
- 4 CE

LCD display 0.01g high precision analytical balance

Feature

- 1.Play with back light LCD
- 2.RS232C communication(optional)
- 3.Toggle between units (g oz ct kg lb)
- 4.Default and auto lock present working mode
- 5.Anti-Radiate, Anti-Static Metal Housing
- 6.Weighing Response time Adjustable
- 7.Counting And Weighing
- 8.Time & Date Setting

Technical Data

| Item No. | Max Cap | Repeat-ability | Linearity | Work Temp | Pan Size | Warm-up Time |
|----------|---------|----------------|-----------|-----------|-----------|--------------|
| ZL-A100 | 100g | ±0.01g | ±0.02g | 5°C-35°C | 128 | 10-20minutes |
| ZL-A200 | 200g | ±0.01g | ±0.02g | 5°C-35°C | 128 | 10-20minutes |
| ZL-A300 | 300g | ±0.01g | ±0.02g | 5°C-35°C | 128 | 10-20minutes |
| ZL-A500 | 500g | ±0.01g | ±0.02g | 5°C-35°C | 128 | 10-20minutes |
| ZL-A600 | 600g | ±0.01g | ±0.02g | 5°C-35°C | 128 | 10-20minutes |
| ZL-A1000 | 1000g | ±0.01g | ±0.02g | 5°C-35°C | 138 | 10-20minutes |
| ZL-A2000 | 2000g | ±0.02g | ±0.04g | 5°C-35°C | 138 | 10-20minutes |
| ZL-A3000 | 3000g | ±0.02g | ±0.04g | 5°C-35°C | 156*156mm | 10-20minutes |

Tel: [+86-769-88055358](tel:+86-769-88055358)

Fax: [+86-769-88055357](tel:+86-769-88055357)

Website: www.zltester.com

E-mail: Zhongli@zltester.com

Address: [Room 722, Dongyuan Mansion, No.159 QiFeng Road, Dongguang city,P.R.C](#)